


International Registry on Cholangiocarcinoma (CCA)



European CCA Registry

Sample Processing Protocol

SAMPLE PROCESSING PROTOCOL

➤ **Serum**

- 1- Collect blood in serum tube
- 2- Keep upright 30 min at Room Temperature
- 3- Centrifuge at 1,500g, 20 minutes, at Room Temperature
- 4- Aliquote the serum in volumes of 1 ml

➤ **Plasma**

- 1- Collect blood in EDTA tube
- 2- Centrifuge at 1,500g, 20 minutes, at Room Temperature
- 3- Aliquote the plasma in volumes of 1 ml
- 3.1- Collect gDNA with **Flexigen DNA Kit** (Automated or Manually)

➤ **Bile**

- 1- Collect bile
- 2- Centrifuge at 1,500g, 20 minutes, at Room Temperature
- 3- Aliquote the bile in volumes of 1 ml

➤ **Urine**

- 1- Collect the urine
- 2- Aliquote 3 ml of urine in volumes of 1 ml
- 2.1- Centrifuge the remainder urine at 3,000g, 5 minutes, at Room Temperature
- 2.2- Aliquote the urine in volumes of 1 ml

➤ **Stools**

- 1- Collect stools
- 2- Freeze at -80°C

➤ **Saliva**

- 1- Collect saliva with **Orogen DNA kit (Genotek)**
- 2- Follow the protocol

➤ **Tumor tissue and surrounding tissue**

▪ **Frozen**

- 1- Cut a piece of non-tumoral and tumoral tissue
- 2- Put them of eppendorf tubes and storage them at -80°C

▪ **Frozen on OCT**

- 1- Cut a piece of non-tumoral and tumoral tissue
- 2- Put each piece in a cryomold with OCT
- 3- Immerse the cryomolds into isopentane at -80°C
- 4- Keep the cryomolds for 5-10 minutes in isopentane
- 5- Leave it drying and storage at -80°C

European CCA Registry (EU-CCA)

- 6- If there is enough tissue, put also a piece of non-tumoral and tumoral tissue in a eppendorf
 - 7- Immerse the eppendorfs into isopentane at -80°C
 - 8- Storage at -80°C
- **Embedded in paraffin**
 - 1- Cut a piece of non-tumoral and tumoral tissue
 - 2- Put each piece in a cassette
 - 3- Immerse the cassettes in formaldehyde 4%
 - 4- Keep in the formaldehyde until the next day
 - 5- Introduce the cassettes in Ethanol 96% during 1h 30 minutes (x2) and then introduce them in Ethanol 96% for 1h
 - 6- Put the cassettes in Ethanol 100% for 1h 30 minutes (x2)
 - 7- Place the cassettes in xilol for 1 h (x3)
 - 8- Soak up the cassettes in paraffin at 55-60°C for 2 h (x2)
 - 9- Put the paraffined cassettes on a “cold plate” to form the block
 - 10- Cut and make the slides